



Delivering Results Through Performance



ADEPT[®]
Hip Resurfacing System




ADEPT[®] Hip Resurfacing System - A Trusted, First Class Solution.




First Class Clinical Data

- Proven low NJR revision and mortality rates
- Over 10 years' data
- Higher function for physiologically young patients
- Early intervention for patients likely to outlive a THR
- Based on 50-year heritage of successful MoM THRs



Optimised Design

- Proven fixation
- Optimised clearance for long-term low wear
- Coverage constant throughout range
- Metallurgy optimised from the 1950s
- Still produced by the original manufacturing specialists



Superior Function

- Constant ROM for all sizes
- Extremely low rate of dislocation compared to THR
- Bone conserving
- Low friction, fluid-film lubrication bearing
- Natural feel and ROM



Simple, Proven Instrumentation

- Optimised sizing and preparation tools
- Simple, reliable neck centring devices
- Simple, reliable cup alignment technique
- Low profile cup introduction
- Safe head cement pressurisation



Implant and Sizing Options

- 2mm head size increments to reproduce natural head diameter
- Minimal cup wall thickness for minimal reaming
- Two cup sizes for every head size
- Extra-fixation cup option

Experience normal again